


# Posture & Torso Alignment (PTA) System Sizing Chart



## Sizing Guidelines



		US – POUNDS / INCHES			METRIC – KG / CM			PTA System
Size Group		Weight	STA	ATK	Weight	STA	ATK	
<b>PEDIATRICS</b>	<b>Infant</b>	15 – 25	8 – 10	7 - 10	7 - 11	20 – 25	18 - 25	q GPTA110
	<b>Pre-School</b>	25 – 45	11 - 13	11 - 12	12 - 20	28 - 33	28 - 31	q GPTA210
	<b>Pediatric</b>	45 – 65	14 - 15	13 - 18	21 - 30	36 – 38	33 - 46	q GPTA310
	<b>Junior</b>	65 – 110	16 - 17	19 - 21	31 - 54	41 - 44	48 - 53	q GPTA410
<b>ADULT</b>	<b>Female Med</b>	110 – 170	18 – 20	22 - 24	55 - 77	46 - 51	56 - 61	q GPTA610
	<b>Female Large</b>	170 – 200	21 - 22	25 - 27	78 - 91	53 – 56	63 - 69	q GPTA620
	<b>Male Med</b>	110 – 170	18 – 19	22 - 24	55 - 77	46 - 49	56 - 61	q GPTA710
	<b>Male Large</b>	170 – 190	20 - 22	25 - 27	78 - 86	50 - 56	63 - 69	q GPTA720

**NEED SIZING HELP? CONTACT CUSTOMER SUPPORT: (970) 239-0344 • CUSTOMERSUPPORT@THERATOGS.COM**